


 Modified Morphological Filtering (MMF) is a filtering technique used to «clean up» the signal with minimal signal distortion


 It allows to obtain reliable results and therefore facilitate the interpretation of the ECG signal




- I. Work on a data set of ECG signals
- 2. Add variable noise to the signals
- 3. Add variable baseline drift to the signals
- 4. Conditioning the signals by MMF algorithm and MF anlgorithm
- 5. Evaluate the performance of the algorithms
- 6. Evaluate the two algorithms by varying the length of the structuring element and the noise applied to the signals



- In order to add noise to the ECG signal, the Additive White Gaussian Noise (AWGN) was chosen for its simplicity and effectiveness in application
- The SNR (Signal to Noise Ratio) is computed randomly to distribute different noises within the data set







- Baseline Drift is simulated by adding a slanted line to a sinusoidal signal
- The function called **GenDrift**, takes care of randomly generating a drift to be applied to all signals in the data set

BASELINE DRIFT APPLIED Drifted ECG signal TO A NOISY SIGNAL Voltage (Hz) Time (sec)




Basic morphological operators:

- **Erosion** (⊕): is a «shrinking» operation. Basically consists in subtracting the structuring element from the input signal, and finding the minimum value of the differences.
- **Dilatation** (\oplus) : is an «expansion» operation. Basically consists in adding the structuring element to the input signal, and finding the maximum value of the sums.


From the first one are derived the following operation:


- **Opening** (o): used to suppress peaks.
- Closing (•): used to suppress pits.





BASELINE DRIFT CORRECTION

- Baseline drift correction is performed by estimating the drift and subtracting it from the signal
- Two structuring elements are used: one for removing peaks (B_o) and the other (B_c) for removing pits left after the previous passage
- B_o and B_c are built as two horizontal line segments of zero amplitude with the lengths that depends on the duration of the characteristic wave and the sample frequency
- Detected badeline drift: $f_b = f_o \circ B_o \bullet B_c$











MMF: NOISE SUPPRESSION

- Noise suppression is performed by processing the signal through a sequence of opening and closing operation concurrently, and then the results are averaged.
- The operations are carried out through a structuring element called B_{pair} , which is composed of two elements, B_1 and B_2 , having the same length to minimize signal distortion
- B_1 has a triangular form used to retain the peaks and valley of the characteristic waves
- B_2 is a line segment used to removing noise


•
$$f = \frac{1}{2} (f_{bc} \bullet B_{pair} + f_{bc} \circ B_{pair}) = \frac{1}{2} (f_{bc} \oplus B_1 \ominus B_2 + f_{bc} \ominus B_1 \oplus B_2)$$








- The performance of the algorithms are evaluated based on the level of noise reduction achieved and on the distortion caused in the signal
- Two main parameters are judged: the Noise-Suppression ratio
 (NSR) and the Signal-Distortion ratio (SDR)
- To give reliable results it is important to minimize the distortion of the ECG signal caused by signal conditioning
- The MMF algorithm retain the significant singular points, therefore is preferred







- Final evaluation performed on the data set as the length of the structuring element grows and the SNR (Signal to Noise Ratio) increases
- For N structuring element (generated through 'GenStrel'), run the MMF and MF algorithm on the entire dataset
- For each mmf-th and mf-th, compute the NSR and the SDR, obtaining N sets one for each structuring element applied, representing the evaluation signal by signal
- For each N-th set calculate the mean, thus obtaining the mean values of the NSR and SDR calculated for the MMF and MF algorithm

